

Your Sorption Science Expert



Surface Measurement Systems

World Leader in Sorption Science





Surface Measurement Systems

World Leader in Sorption Science



 Manufacturing  Lab Facilities  Sales & Support  Regional Service Support

Instruments – Product lines



DVS

IGC SEA



Surface Measurement Systems
World Leader in Sorption Science

New Products



BTA



iGC-SEA
Inverse Gas Chromatography-
Surface Energy Analyzer



DVS
Carbon


BTA Frontier

Specifications



- **Three different gas or vapor inlets/sources:** CO₂, H₂O and VOCs
- Generation of water vapor from 0-90%RH
- **Small amounts of sample needed (~10 to 1000 mg)**
- Vapor/gas sensors are located before & after sample column
 - CO₂, Water vapor, PID sensor- organic species: 1ppb to 100ppm and Thermal conductivity detector (TCD)
- **In-situ sample conditioning/regeneration up to 300 °C**

Breakthrough curve: crucial parameters




Measure true single/multicomponent sorption by determining uptakes of each.

Obtain crucial process parameters like breakthrough and equilibration time.

Investigate kinetics, diffusion, flow rate effects, etc.

Case study: Impact of H₂O on CO₂ adsorption



Influence of water on the CO₂ uptake: 20% RH and 3.5% vol CO₂ on **Zeolite 13X**. Three sequential breakthrough experiments (activation not shown).

- Explore the **competitive effect of water on the CO₂ adsorption** of porous materials
- Characteristic “roll-up” effect is observed, as the water front is replacing adsorbed CO₂
- **Near-perfect repeatability** over three cycles

High Temperature IGC-SEA



New high temperature oven: **30 °C - 500 °C**


Conduct surface energy analysis under industrial conditions

Fully-automated iGC system

Patented headspace injection system with humidity generator

Purpose-built control & analysis software with CFR capabilities

Dispersive Surface Energy of Zeolite 13x



At 300° & 350 °C is homogeneous, the surface energetically uniform.

At 400 °C is shows some degree of heterogeneity.

DVS Carbon



Main Features

- CO₂ gravimetric sorption analyzer
- Water vapor sorption analyzer
- Co-adsorption of CO₂ and water vapor
- Easy access to kinetics, sorption/desorption and hysteresis measurements
- Patented Speed of Sound for CO₂ measurement
- True0 drying at 0.0% RH
- Dual configuration for broad range of CO₂ concentrations*
 - High (0 -98%)
 - Low (e.g 400 ppm)
 - *Requires a cylinder of CO₂ to be connected to the DVS instrument

DVS

Thank you for having us, the opportunity
and see you in 2024!



Surface Measurement Systems
World Leader in Sorption Science



DVS Instruments*

Water Vapour Only



DVS
Intrinsic Plus

Water vapor
20-40 °C
Small footprint



DVS
ADVENTURE

Water vapour
5-85 °C

Water and Organic Vapours



DVS
RESOLUTION

Water and organic vapour
5 – 85 °C



DVS
DISCOVERY

Water or organic vapour
10 – 70 °C (Dual balance)



DVS
ENDEAVOUR

Water and organic vapor
10 – 70 °C (Five balances)

*Flow based - atmospheric pressure systems